Silentflow 2

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4.2 Procedures and Disinfecting Procedures:

The following procedures are allowed:

Cleaning by hand

Remove all residue with warm water and mild detergent, just as you would wash your dishes. Rinse all parts carefully with clean water.

Dishwasher

Wash cycle at maximum 65° C. Afterwards, rinse all parts carefully with clean water.

Disinfection*

Clean all parts as described in "Cleaning by hand" and then dip the parts in a diluted disinfectant solution. After disinfection, rinse all parts thoroughly with **distilled** water.

Dry all parts thoroughly after cleaning. Check the components for any remaining dirt or damage. Repeat cleaning if necessary. If damaged, Silentflow 2 must be replaced. Assemble Silentflow 2 as described above.

* We recommend **Gigasept FF**. Please be sure to observe the manufacturer's instructions. Tolerance to this disinfection agent has been proven. We recommend that you wear suitable gloves when disinfecting the equipment (e.g. household or disposable gloves).

4.3 Notes on disinfection

It may be necessary to disinfect the parts after each use by patients with a weakened immune system or a certain clinical background. Please consult your physician in this case.

4.4 Disinfection with change of patient

If the exhalation system is to be used for another patient, it must be hygienically treated beforehand. For this purpose, dismantle and disinfect all parts as described above or use a new exhalation system.

5 Maintenance

Silentflow 2 is maintenance-free.

When Silentflow 2 can no longer be used, dispose of it with your regular hospital waste disposal or with normal household waste removal.

6 Supply Schedule

Silentflow 2	WM 23600
consists of:	
– Silentflow 2	
– Instructions for use DE, EN, FR	WM 16656
Other available products:	
Instructions for use IT, NL, ES	WM 16776
Instructions for use DA, NO, SV	WM 16815


All parts listed are available as replacement parts.

7 Technical data

Device class according to 93/42/EEC:	Ila
Connection to mask:	Cone 22 mm diam. (EN ISO 5356-1)
Connection to hose: suitable for hoses made of	Ø 21.5 mm, 1° taper, elastomer with an interior connection diameter of 20 to 21 mm
Therapy pressure range:	4 hPa to 40 hPa
Weight:	15 grams
Temperature range:	
– Operating:	+5°C to + 40°C
– Storage:	-20°C to + 70°C
Flow resistance at 50 l/min	0,04 hPa
Flow resistance at 100 l/min	0,18 hPa
Quoted two-number noise emission value to ISO 4871	
– sound pressure level	16 dB(A)
– sound power level	24 dB(A)
(uncertainty factor: 3 dB(A))	
Usage life:	6 months, with careful maintenance up to 12 months, depending on the frequency of use and cleaning
Material:	PC (polycarbonate)

CE 0197

Pressure-flow characteristic curve
In the pressure-flow characteristic curve, the outlet flow is shown in dependence with therapy pressure.

8 Warranty

Weinmann guarantees that the product, when used in accordance with requirements, will remain free from defects for a period of 6 months from date of purchase.

Claims against the warranty can be made only when accompanied by the sales receipt, which must show salesperson and date of purchase.

We offer no warranty in the case of:

- Disregard of usage instructions
- Operating errors
- Improper use or handling
- Third-party intervention by non-authorized persons for the purpose of device repair
- Acts of God
- Transport damage as a result of improper packaging of returned items
- Operational and normal wear and tear
- Failure to use original spare parts.

Weinmann is not liable for consequential harm caused by a defect if it is not based on intention or gross negligence. Weinmann is also not liable for minor physical injury to life or limb resulting from negligence.

Weinmann reserves the right to decide whether to eliminate defects, to deliver a defect-free item or to reduce the purchase price by a reasonable amount.

If Weinmann rejects a claim against the warranty, it assumes no expense for transport between customer and manufacturer.

Implied rights to claim under warranty remain unaffected by these changes.

Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG declares herewith that the product complies fully with the respective regulations of the Medical Device Directive 93/42/EEC. The unabridged text of the Declaration of Conformity can be found on our website at: www.weinmann.de

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1 Device description

Legend

- 1 Connecting sleeve
- 2 Middle section
- 3 Opening
- 4 Connection to hose
- 5 Connection to mask

Purpose of Use

The exhalation system Silentflow 2 is an accessory that is used with nasal masks in the treatment of sleep apnea and in non-invasive ventilation of patients with respiratory insufficiency.

Silentflow 2 makes it possible for you to exhale quietly during therapy. Silentflow2 prevents the exhaled air from being inhaled again.

Functional description

Silentflow 2 connects the breathing hose of the therapy device to the nasal mask.

The inhaled air flows through the mask-hose system. The exhaled air escapes through an opening between the middle section and connecting sleeve.

2 Safety Requirements

Special markings on the device

- Storage temperature

–20 °C

+70 °C

- Outlet for exhaled air

Safety Instructions

For your own protection and that of your patient and according to the requirements of the European Union Guidelines 93/42/EEC, please observe the following:

- Please read these instructions for use carefully. They are a part of the device and have to be available at all times.
- Use only nasal masks that do not have their own exhalation openings.
- Do not use masks that cover the entire nose and mouth.
- Silentflow 2 has to be directly connected to the mask.
- Silentflow2 is set up with a therapy pressure range between 4 and 40 hPa. Silentflow2 cannot be used with lower or higher therapy pressures.
- Please note that patients with limited spontaneous breathing must be monitored at all times and that an alarm for negative pressure and/or leaks has to be activated on the ventilation device.
- Patients with obstructive and restrictive pulmonary illnesses should be individually examined, for example with regard to CO₂ rebreathing, during usage of Silentflow.
- If third-party items are used, functional failures may occur and fitness for use may be restricted. Biocompatibility requirements may also not be met. Please note that in such cases, any claim under warranty and liability will be voided if neither the accessories nor genuine replacement parts recommended in the instructions for use are used.
- Make sure that the connections between Silentflow 2 and the mask and between Silentflow 2 and the breathing hose are firm and air-tight. The connections must not be allowed to loosen during therapy.
- Do not close any openings on Silentflow 2.
- Place the nasal mask on your face for any length of time only when the therapy device has been switched on. One of the purposes of the therapy device air flow is to remove the exhaled air from the mask and hose. When the therapy device is not switched on, there is a risk of rebreathing the exhaled air if you breathe through your nose.
- Do not use any electrically conductive breathing hoses.
- Clean the exhalation system as described in Point 4 "Hygienic preparation". Unsuitable cleaning methods, unsuitable detergents and insufficient cleaning frequency can lead to infection or bacterial contamination and impair the function of the exhalation system.

3 Operation

Putting into operation

Connect Silentflow 2 to the breathing hose. Ensure that the air outlet is directed toward the hose.

Insert the Silentflow 2 into the tapered connection of the mask. You can make your therapy more comfortable by observing the following points:

Six arrows indicate the area on the connecting hose where exhaled air escapes. There is a segment in between, in which the exhalation opening is closed. Insert the exhalation system into the tapered connection of the mask so that this segment is directed toward your body. This will prevent an uncomfortable draft onto your body.

After use

Disconnect Silentflow 2 from the mask and breathing hose. Clean Silentflow 2 as described in Point 4 "Hygienic preparation".

Separating the exhalation system

Unlock the middle section and connection hose by applying slight pressure to the hook and by pulling the pieces apart.

Assembling the exhalation system

Fit the parts together. Ensure that the lock engages (clicks into place). When Silentflow 2 is correctly assembled, the middle section and hose can be rotated in opposite directions from one another.

Combination with other devices

This exhalation system can be used with all devices for treatment described under "Purpose of Use" if the following conditions are fulfilled:

- When Silentflow 2 is combined with other devices, ensure that the therapy pressure and the pressure constancy in the nasal mask are as prescribed by your physician. Ask your physician or dealer whether the correct therapy pressure will be achieved in the mask using the desired equipment combination. Have your physician or dealer check the settings of your therapy device.
- Follow the directions in the user's manuals of the devices you want to combine with Silentflow 2.
- Use a nasal mask with a tapered connection in accordance with EN ISO 5356-1, e.g., SOMNOmask.
- Use a breathing hose made of elastomer with a connection diameter of 20 to 21 mm.

Use of heated and unheated humidifiers

Silentflow 2 can be used with a heated or an unheated humidifier. Ensure that Silentflow 2 is dry before starting treatment. Never select a degree of humidification which is so high that water forms in the breathing hose or exhalation system. Refer to the user's manuals for your humidifier and therapy device.

4 Hygienic preparation

4.1 Intervals

A dirty exhalation system can make exhalation difficult and prevent gas exchange. Therefore, it is important to clean Silentflow 2 to assure trouble-free functioning of the mask.

Clean Silentflow 2 after each use. To do this, dismantle the exhalation system as described above.

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