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GEM is an electrical fed compressor with aerosol therapy atomiser for Home-Care.
The instrument is designed for easy transport and handling and is recommended for atomising antibiotics and bronchodilator drugs. The high thermal insulated plastic body complies with European Safety Rules. The medical device is designed for continuous use. Device designed for non-continuous use (Ton: 20 min / Toff: 40min).

GENERAL WARNING

READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

DRUG ADMINISTRATION MUST BE UNDER MEDICAL CONTROL

THE INSTRUMENT MUST NOT BE DISASSEMBLED, FOR A TECHNICAL SERVICE ALWAYS CONTACT
MOBIAK S.A

REMOVE THE DEVICE FROM CARRYING BAG BEFORE USE

IMPORTANT SAFETY RULES

1. On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and / or peeling of the power supply cable. In these cases don't connect the plug to the electric socket. Carry out these controls before each use;
2. before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. Never leave the appliance inserted if not necessary disconnect the plug from the mains power supply when it is not being used;
4. Respect the safety regulations indicated for electrical appliances and particularly:
 - Only use original accessories and components provided by the manufacturer MOBIAK S.A to guarantee the highest efficiency and safety of the device;
 - Never immerge the appliance into water;
 - Position the appliance on flat stable surfaces;
 - Position the device in a way that the air inlets on the back aren't obstructed;
 - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
 - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
 - The use of this device by children and / or incompetent person always requires the careful surveillance of an adult in possession of their full mental faculties;
 - The medical device, and most of all the nebulae, must be kept out of children's reach as it contains small parts that could be swallowed;
 - Don't leave the appliance connected to the power supply socket when not in use;
 - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
 - Store and use the device in places protected against the weather and far from any sources of heat. After each use, it is recommended to store the device in its own box away from dust and sunlight.
 - In general, it is inadvisable to use single or multiple adapters and/or extensions. Should their use be necessary, you must use ones that are in compliance with safety regulations, however, taking care not to exceed the maximum power supply tolerated, which is indicated on the adapters and extensions.
5. For repairs, exclusively contact MOBIAK S.A technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device;
6. This medical device must be destined exclusively for the use for which it has been designed as described in this manual. It must therefore be used as an aerosol therapy system. Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force;
7. The medical device requires special precautions regarding electromagnetic compatibility and must be installed and used in accordance with the information provided with the accompanying documents: the GEM device must be installed and used away from mobile and portable RF communication devices (mobile phones, transceivers, etc.) that may interfere with the said device.
8. Store the accessories out of reach of children. Children and people with learning difficulties must only use the medical device under the strict supervision of an adult with full mental faculties. Keep the ampoule out of reach of children under 36 months as it contains small parts that may be swallowed accidentally. Never leave the device unattended in places accessible to minors and / or the disabled

9. **WARNING:** None of electric or mechanical parts have been designed to be repaired by customers or end-users. Don't open the device, do not mishandle the electric / mechanical parts. Always contact MOBIAK S.A technical assistance.
10. Do not leave the device unattended in places accessible to children and / or persons not in full possession of their mental faculties as there is a risk of strangulation with the air tube;
11. The medical device may come into contact with the patient via the nebuliser / masks / mouthpiece and / or nosepiece, components compliant with the requirements of regulation ISO 10993-1: therefore, no allergic reaction and skin irritation may occur.
12. The product and its parts are biocompatible in accordance with the requirements of regulation EN 60601-1.
13. Operation of the device is very simple and therefore no further explanations are required other than those indicated in the following user manual.
14. Using the device in environmental conditions different than those indicated in this manual may harm seriously the safety and the technical characteristics of the same;
15. The materials used to contain the drugs are made with highly stable thermoplastic polymers that are resistant against chemicals. Such materials were tested with commonly used drugs (Salbutamol, Beclomethasone dipropionate, Acetylcysteine, Budesonide, Ambroxol) and no interaction phenomenon was observed. Interactions cannot however be excluded given the variety and the continuous evolution of the drugs that are used. Remember to:
 - To consume the drugs as quickly as possible after opening its package;
 - To avoid keeping the drug in the tray-like container for too long and to clean it immediately after every application;
 - If the tray-like container presents any abnormal situation (such as softening or cracks), do not introduce any solution and do not proceed with the inhalation. Contact the technical service and describe the methods and type of drugs used.
16. Remember to:
 - Only use this device with medicines prescribed by your doctor;
 - Carry out the treatment only using the accessory indicated by the doctor according to the pathology.

Under certain failure conditions the temperature of the casing may become hot and there may be a risk of burns if you touch those parts. In any case, the temperatures do not exceed the limit of 105°C (ref. Interpretation Sheet IEC 60 601 -1).

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

At the end of its working life, the product must not be disposed of as urban waste. It must be taken a special local authority differentiated waste collection or to a dealer providing this service. Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant saving in energy and resources. As a reminder of the need to dispose of separately the product is marked with a crossed-out wheeled dustbin.

RULES FOR RETURNING AND REPAIRING

UNDER NEW EUROPEAN RULES, MOBIAK S.A. REQUIRES THE FOLLOWING PROCEDURES TO BE CARRIED OUT TO PROTECT THE INSTRUMENT AND THE SAFETY OF ALL WHO COME IN CONTACT WITH IT.

MOBIAK S.A. warrants its product for 60 months after purchasing date.

Before returning an instrument for repair, the external surfaces and all accessories **MUST** be carefully disinfected with a cloth soaked in methylated spirits or hypochlorite-based solution. The instrument and accessories should then be placed in a bag with a note outlining the disinfection undertaken.

Failure to follow this procedure will result in the instrument being returned to the purchaser unrepaid.

Instruments returned for repair **MUST** be accompanied by a description of the problem. MOBIAK S.A. will not be responsible for damage caused through improper use. To avoid such damage, please read the instruction carefully.

Where MOBIAK S.A. determines that an instrument is faulty, a replacement will be provided only if a **SALES RECEIPT** and **STAMPED GUARANTEE** are provided. MOBIAK S.A. will not be responsible for damage accessories. These may be replaced at the customer's expense.

SIMBOLOGY

	Class 2 isolation equipment
CE0653	CE marking in conformity with EC directive 93/42/CEE and subsequent changes
	General warnings and/or specifications
	Consult the instruction manual
	Keep in a cool, dry place
	Conservation temperature: - 25 ÷ 70° C
	Atmospheric Pressure
	Applied part type BF (Nebulizer, mouthpiece, nosepiece, pediatric mask and Adult mask)
~	Alternate Current
Hz	Mains frequency
I	ON
O	OFF
LOT	Batch production
SN	Serial number
REF	Model / Ref number
	MOBIAS S.A 96-98 M. Mpotsari STR. Chania, Crete, Greece

The manufacturer cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its components be damaged due to accident or misuse.


Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC (and subsequent changes) and its normatives.

The accessories that accompany the device are compatible with EN 13544-1.

The device cannot be used in an anesthesia circuit or respirator circuit.

TECHNICAL CHARACTERISTICS

TIPOLOGY (MDD 93/42/EEC)	Class IIa Medical Device	
MODEL	GEM	
POWER FEEDING	230V~ / 50Hz	
POWER CONSUMPTION	170 VA	
MAX PRESSURE	250 kPa (2.5 Bar)	
MAX AIR FLOW	14 l/min	
OPERATING PRESSURE	110 kPa (1.10 Bar)	
OPERATING AIR FLOW	4.8 l/min at 110 kPa	
NEB-RATE	0.35 ml/min (with 4ml of 0.9% NaCl solution)	
MMAD	2.44 μ m	
GSD	2.87	
PESO	1.45 Kg	
DIMENSIONS	130 (L) x 89 (H) x 146 (P) mm	
NOISE LEVEL	60 dB (A)	
DUTY CYCLE	Ton: 20min / Toff: 40min	
MIN CAPACITY NEBULIZER	2ml	
MAX CAPACITY NEBULIZER	6ml	
WORKING CONDITION	Room Temperature:	5 $\div$ 35°C
	Room humidity percentage	15 $\div$ 93% RH
	Atmospheric pressure:	800 $\div$ 1060 hPa
CONSERVATION AND TRANSPORT CONDITIONS	Room Temperature:	-25 $\div$ 70° C
	Room humidity Percentage:	0 $\div$ 93% RH
	Atmospheric pressure:	500 $\div$ 1060 hPa

MMAD = Mass Median Aerodynamic Diameter

GSD = Geometric Standard Deviation

NB: The measures and curves are not valid for the high viscosity suspension drug.

MAINTENANCE

The GEM aerosol does not need maintenance or lubrication.

With regard to training, given the information contained in the user manual and since it is easy to understand the said device, it doesn't appear to be necessary. It is necessary to check functioning and instrument before every use. Unpack the instrument and always check integrity of plastic parts and feeding cable., they might have been damaged during previous use.

Connect cable to electrical network and turn switch on. Close the compressor mouthpiece with a finger and verify that loud noises are not present, these can indicate wrong functioning. Verify that the atomiser is not damaged by previous use (it was badly put away or badly knocked).

A protection fuse (F 2A L 250V) not reachable from exterior protects the instrument.

For fuse replacing, please make reference to manufacturer technical personnel.

MOBIAK S.A will provide upon request electric diagrams, components list, description, setting instructions and any other information that can help the technical assistance staff for product repair.

Fault type	Cause	Solution
1. Low Nebulization	Clogged Nebulizer Tank	Clean and disinfect the nebulizer tank as explained in the instruction manual
2. Low Nebulization	Clogged Nebulizer Tank	If cleaning was not succesful change cruet
3. Absence of Nebulization	Clogged Nebulizer Tank	Check that the nebulizer contains medication; Make sure that the nebulizer is not clogged; Check the connection between the compressor air outlet port and the accessories
4. Slow Nebulization	Highly dense drug	Dilute drug in physiological liquid
5. The device doesn't work	Power Cable is damaged External power source failure	Replace the power cable Check the external power source
6. Noisy Device	Extended use	Call retainer or manufacturer MOBIAK S.A
Fault 1 - 2 - 3 - 4 - 5 - 6	No solution with previous items	Call retainer or manufacturer MOBIAK S.A

If the unit doesn't nebulizer once the above conditions have been checked, we suggest to contact your dealer or technical service MOBIAK S.A.

BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT CA TECHNICAL SERVICE. MOBIAK S.A DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.

ANY MINIMAL MODIFICATION/REPAIR ON THE DEVICE VOIDS THE WARRANTY AND DOES NOT GUARANTEE THE COMPLIANCE WITH THE TECHNICAL REQUIREMENTS PROVIDED BY THE MDD 93/42/EEC DIRECTIVE (AND SUBSEQUENT CHANGES) AND ITS NORMATIVES

CLEANING DEVICE

Use a soft dry cloth with not – abrasive and not – solvent detergents. The device's plug must be removed from the wall socket before proceeding with any cleaning procedures.

PARTICULAR CARE SHOULD BE TAKEN TO ENSURE THAT THE INTERNAL PARTS OF THE EQUIPMENT DO NOT GET IN TOUCH WITH LIQUIDS. NEVER CLEAN THE EQUIPMENT WITH WATER.

CLEANING OF ACCESSORIES

Before using the device, the manufacturer advises you to clean the accessories.

Proceed as follows to disinfect the accessories:

1. Turn the upper part of the nebulizer in an anti-clockwise direction;
2. Disconnect the internal pipette at the base of the nebulizer using the fingers;
3. After using the appliance, disassemble the nebulizer and clean all parts in warm water, rinse carefully and remove excess water using a soft cloth and leave to dry in a clean place.
4. To disinfect the accessories (except the tube) use denatured alcohol or hypochlorite solutions easily obtainable from chemists. Before using after cleaning, rinse the accessories in warm water until all traces of the disinfectant have been eliminated, dry and store in a dry and dust free place.


DO NOT BOIL OR PUT IN AUTOCLAVE MASK AND AIR TUBE

After each treatment clean thoroughly each component of the nebulizer (except air tube) removing medication residual and possible impurities. Clean all parts in warm water. Rinse thoroughly making sure that all deposits are washed away and let dry.

STANDARD ACCESSORIES

ACCESSORIES
HI-FLO KIT - REF 0148006/7 (Nebulizer HI-FLO, Adult Mask, Pediatric Mask, Air Tube and Mouth-piece, Nosepiece)
Air Filter (+ n°2 spare parts)
Carrying Bag

For each individual patient it's recommended to use the nebulizer for 6 months or for a maximum of 120 treatments. The nebulizer must be replaced after a long period of inactivity, if it is deformed or broken, or if the nebulizer nozzle is blocked by dry medicine, dust, etc. Only use the original nebulizer supplied by MOBIAS S.A with the device. Use the "nose piece" accessory only if expressly indicated by your doctor and paying attention NEVER to introduce inside the nose the nasal bifurcation, but only bring it as close as possible.

In the presence of infection or microbial contamination prone pathologies, we recommend using your personal accessories and nebulizer (always consult your doctor).

The device is equipped with a filter that removes any impurities from the air that was sucked in by the compressor. The air filter must be replaced every 25 hours of functioning or when it result particularly worn. For replacement, lift the filter and replace with a new one. Only use original MOBIAS S.A filter.

The mask and tube must be replaced as soon as the materials they are made of show signs of deterioration.

Service life of the device: More than 1500 hours of operation (or 5 years) in accordance with the standard conditions of testing and operation. Shelf life: maximum 5 years from the date of manufacture.

INSTRUCTION FOR USE

The device must be checked before each use in order to detect malfunctions and / or damage caused by transport and / or storage. During the inhalation must sit in an upright and relaxed position at a table and not in an armchair, to avoid compressing the airways and therefore compromising the effectiveness of the treatment. It is recommended not to keep the device in your hands and / or to avoid prolonged contact with the body of apparatus.

WARNING: Put the device on a flat and stable surface in order not to block the cooling vents on the sides of the device.

1. Extract the power supply cable and insert the plug into the mains socie. It is recommended to unwind the entire length of the power supply cable to prevent dangerous overheating. If the power supply cable is damage and must be replaced contact MOBIAM S.A technical service.
2. Open the nebulizer by unscrewing the lid;
3. Pour the medicine prescribed by the doctor into the nebulizer;
4. Re-close the nebulizer, re-screwing the lid;
5. Connect air pipe to rhe air exit well;
6. Connect the other end of the pipe to the connection in the lower part of the nebulizer;
7. Connect the desidered accessory to the nebulizer (pediatric mask, adult mask, mouthpiece, nosepiece);
8. Ensure that the supplied air filter is present;
9. Press switch on position I to proceed with nebulization;
10. On completing of nebulization, press the switch on position 0 and remove the plug from the socket;
11. Wash the nebulizer and its accessories as indicated in the cleaning charter;
12. Place the cable and accessories into the box.

Always use the nebulizer facing upwards so that substances and / or medicines cannot escape from the nebulizer during the normal use.

WARNING: The power supply cable plug is the element of separation from the electrical mains system: even if the units equipped with a special on / off switch button, the power supply plug must be kept accessible once the device is in use so as to allow a further method of disconnection from the mains supply system.

NEVER INHALE IN HORIZONTAL POSITION

NEVER BENS THE NEBULIZER OVER 60°

RISK OF ELECTROMAGNETIC INTERFERENCE AND POSSIBLE REMEDIES

This section contains information regarding the conformity of the compliance with the EN 60601-1-2 Standard. The GEM is an electro - medical device that requires particular precautions regarding electro -magnetic compatibility and which must be installed and commissioned according to the electro-magnetic compatibility information supplied. Portable and mobile radio communication devices (mobile phones, transceivers, etc.) may interfere with the medical device and should not be used in close proximity with, adjacent to or on top of the medical device. If such use is necessary and unavoidable, special precautions should be taken so that the electro-medical device functions properly in its intended operating configuration (for example, constantly and visually checking for the absence of anomalies or malfunctions). The use of accessories, transducers and cables different to those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as spare parts, can lead to an increase in emissions or in a decrease of the immunity of the device or system. The following tables supply information regarding the EMC (Electromagnetic Compatibility) characteristics of the electro-medical device.

Guidance and manufacturer's declaration – Electromagnetic Emissions		
The aerosol GEM is intended for use in the electromagnetic environment specified below. The customers or the user of the aerosol GEM should assure that it's used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The aerosol GEM only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The aerosol GEM can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes or environments feeds to you from batteries.
Harmonic emissions EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – Immunity Emissions			
The aerosol GEM is intended for use in the electromagnetic environment specified below. The customers or the user of the aerosol GEM should assure that it's used in such an environment.			
Immunity Test	Level indicated by the EN 60601 -1-2	Compliance Level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 8kV on contact ± 15kV in air	The device doesn't change its state	Floors should be wood, concret or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst EN 61000-4-4	± 2kV power supply lines ± 1kV for input / output lines	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Surge EN 61000-4-5	±0,5 and ±1kV differential mode	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Loss of voltage, brief voltage interruptions and variations EN 61000-4-11	5%U _T (>95% dip U _T) for 0.5 cycle 40%U _T (>60% dip U _T) for 5 cycle 70%U _T (>30% dip U _T) for 25 cycle <5%U _T (>95% dip U _T) for 5 sec	-	Mains power quality should be that of a typical commercial environment or hospital If the user of the aerosol GEM request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field EN 61000-4-8	30A/m	The device doesn't change its state	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Nota U _T is the value of the power supply voltage			

Guidance and manufacturer's declaration – Immunity Emissions			
The aerosol GEM is intended for use in the electromagnetic environment specified below. The customers or the user of the aerosol GEM should assure that it's used in such an environment.			
Immunity Test	Level indicated by the EN 60601 -1-2	Compliance Level	Electromagnetic environment - guidance
Conducted Immunity EN 61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting devices)	$V_1 = 3 \text{ V rms}$	The portable and mobile RF communication devices, including cables, must not be used closer to the GEM device, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance $d = [3.5 / V_1] \sqrt{P}$ $d = [12 / E_1] \sqrt{P} \text{ from 80 MHz to 800MHz}$ $d = [23 / E_1] \sqrt{P} \text{ from 800 MHz to 2.5 GHz}$ Where P is the maximum nominal output voltage of the transmitter in Watt (W) depending on the manufacturer of the transmitter and the recommended separation distance in metres (m). The intensity of the field from the fixed RF transmitters, as determined by an electro-magnetic study of the site^{a)}, could be lower than the level of conformity of each frequency interval ^{b)}. It is possible to check for interference in proximity to devices identified by the following symbol: 
Radiated Immunity EN 61000-4-3	10V/m 80MHz to 2.7GHz (for non life-supporting devices)	$E_1 = 10 \text{ V / m}$	

Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied

Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by reflection from buildings, objects and people.

a) The field intensity for fixed transmitters such as the base stations for radiotelephones (mobile and cordless) and terrestrial mobile radio, amateur radio devices, radio AM and FM transmitters and TV transmitters can not be theoretically and accurately foreseen. To establish an electro-magnetic environment generated by fixed RF transmitters, an electro-magnetic study of the site should be considered. If the field intensity measured in the place where the device will be used surpasses the above mentioned applicable level of conformity, the normal functioning of the device should be monitored. If abnormal performance arises, additional measures such as changing the device's direction or positioning may be necessary.

b) The field intensity on an interval frequency of 150 kHz to 80 MHz should be less than 10 V/m.

Recommended separation distance between portable and mobile radio -communication devices and the monitor			
The aerosol GEM surgical aspirator is intended to operate in an electro-magnetic environment where RF irradiated interferences are under control. The client or operator of the GEM device can help prevent electro-magnetic interference by keeping a minimum distance between the portable and mobile RF communication devices (transmitters) and the GEM device, as recommended below, in relation to the radio-communication maximum output power.			
Maximum nominal output power of the Transmitter W	Separation distance from the frequency transmitter (m)		
	150 kHz to 80 MHz $d = [3.5 / V_1] \sqrt{P}$	80 MHz to 800 MHz $d = [12 / E_1] \sqrt{P}$	800 MHz to 2.5 GHz $d = [23 / E_1] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters with a maximum nominal output power not shown above, the recommended separation distance in metres (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum nominal output power of the transmitter in Watt (W) depending on the transmitter's manufacturer.			
Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied			
Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by the reflection from buildings, objects and people.			

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Notes

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Notes

[illegible]

Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There are no margins, text, or other markings on the paper.