

Κατάλογος συμβατότητας για μάσκες/συσσκευές

Ελληνικά

Το παρόν έγγραφο παραθέτει τις μάσκες που είναι συμβατές με τις συσκευές θετικής πίεσης αεραγωγών της ResMed και τη ρύθμιση της μάσκας, εφόσον απαιτείται. Για πλήρεις λεπτομέρειες σχετικά με τη σωστή χρήση αυτών των συσκευών, παρακαλείστε να ανατρέξετε στις αντίστοιχες οδηγίες χρήσης.

Σημείωση: Δεν είναι όλα τα προϊόντα διαθέσιμα σε όλες τις γεωγραφικές περιφέρειες.

Μάσκες για ολόκληρο το πρόσωπο 1

Συσκευή	AcuCare F1-4	AirFit™ F40 ⁴	AirFit™ F30i ¹ / AirFit™ F30i Non Magnetic ¹	AirFit™ F30	AirTouch™ F20 / AirTouch™ F20 for Her	AirFit™ F20 / AirFit™ F20 for Her / AirFit™ F20 Non Magnetic	AirFit™ F10 / AirFit™ F10 for Her	Mirage™
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	X	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	X/E
AirCurve™ 11 ASV PaceWave	X	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	X/E
AirSense™ 11	X/E	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	X/E
AirMini™	X	X	X	✓	✓	✓	X	X
AirStart™ 10 / AirSense™ 10	X/E	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	X/E
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	X/E	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	X/E
AirCurve™ 10 CS Pacewave	X/E	Pillows	Full Face	Full Face	Full Face	Full Face	Full Face	X/E
AirCurve™ 10 CS-A Pacewave	X/E	X/E	X/E	Full Face	Full Face	Full Face	Full Face	X/E
Lumis™	X/E	Pillows ³	Full Face ³	Full Face ³	Full Face	Full Face	Full Face	X/E
S9 AutoSet™ / S9 Elite™	X/E	Pillows	Full Face	X/E	Full Face	Full Face	Full Face	X/E
S9 Escape™	X/E	Pillows	Full Face	X/E	Full Face	Full Face	Full Face	X/E
AutoSet CS™ 2	X/E	X/E	X/E	X/E	X/E	X/E	FULL FACE	FULL FACE
AutoSet CS™ (S9)	X/E	Pillows	X/E	X/E	Full Face	Full Face	Full Face	X/E
AutoSet CS™ (S9 PaceWave)	X/E	Pillows	X/E	X/E	Full Face	Full Face	Full Face	X/E
AutoSet CS-A™ (S9 PaceWave)	X/E	X/E	X/E	X/E	Full Face	Full Face	Full Face	X/E
C-Series Tango™	X/E	X/E	X/E	X/E	X/E	X/E		Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.
S7™ Lightweight	X/E	X/E	X/E	X/E	X/E	X/E		Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.
S7 Elite™ / S7 AutoSet Spirit	X/E	X/E	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL
S7 AutoSet Respond™	X/E	X/E	X/E	X/E	X/E	X/E	MIRAGE	STANDARD
S8™ Auto 25	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Spirit™	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Vantage™	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 Escape™	X/E	X/E	X/E	X/E	X/E	X/E		Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H ₂ O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.

Συσκευή	AcuCare F1-4	AirFit™ F40 ⁴	AirFit™ F30i ¹ / AirFit™ F30i Non Magnetic ¹	AirFit™ F30	AirTouch™ F20 / AirTouch™ F20 for Her	AirFit™ F20 / AirFit™ F20 for Her / AirFit™ F20 Non Magnetic	AirFit™ F10 / AirFit™ F10 for Her	Mirage™
S8 Escape™ II	X/E	X/E	X/E	X/E	Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.			
S8 Elite™ / S8 Elite™ II	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Spirit™ II	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 Compact™	X/E	X/E	X/E	X/E	X/E	X/E	Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.	
VPAP™ III / VPAP™ III ST	X/E	X/E	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL
VPAP™ III ST-A	X/E	X/E	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL
VPAP™ IV / VPAP™ IV ST	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ Adapt SV	X/E	X/E	X/E	X/E	X/E	X/E	FULL FACE	FULL FACE
VPAP™ Adapt (S9)	X/E	Pillows	Full Face	X/E	Full Face	Full Face	Full Face	X/E
VPAP™ Adapt (S9 PaceWave)	X/E	Pillows	Full Face	X/E	Full Face	Full Face	Full Face	X/E
VPAP™ Malibu	X/E	X/E	X/E	X/E	X/E	X/E	FULL FACE	FULL FACE
VPAP™ Auto 25	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
Auto 25 (S9)	X/E	Pillows	Full Face	X/E	Full Face	Full Face	Full Face	X/E
VPAP™ S / VPAP™ ST	X/E	X/E	X/E	X/E	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ S (S9) / VPAP™ ST (S9)	X/E	Pillows	Full Face	X/E	Full Face	Full Face	Full Face	X/E
VPAP™ ST-A (S9)	X/E	X/E	X/E	X/E	Full Face	Full Face	Full Face	X/E
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	X/E	X/E	X/E	X/E	Full Face	Full Face	Full Face	X/E
S9 VPAP™ Tx	X/E	Pillows ²	Full Face ²	X/E	Full Face	Full Face	Full Face	X/E
VPAP™ Tx	X/E	X/E	X/E	X/E	X/E	X/E	FULL FACE	X/E
Stellar™	X/E	X/E	X/E	X/E	FULL FACE	FULL FACE	FULL FACE	X/E
VS Ultra™ / VS Integra™	✓	X/E	X/E	X/E	X/E	X/E	✓	X/E
VS III™	X/E	X/E	X/E	X/E	X/E	X/E	✓	X/E
Astral™	FULL FACE	X/E	X/E	X/E	FULL FACE	FULL FACE	X/E	X/E

¹ Το SmartStop ενδέχεται να μην λειτουργεί αποτελεσματικά όταν αυτή η μάσκα χρησιμοποιείται με ορισμένες συσκευές CPAP ή δύο επιπέδων.

² Εξαιρούνται οι καταστάσεις λειτουργίας PAC και iVAPS. Συμβατό μόνο για το εύρος πίεσης 4 – 25 cm H2O.

³ This combination of mask and device has been determined safe for use based on the evaluation of accuracy of monitoring of pressure, volume, minute ventilation and respiratory rate, and its alarms if applicable. Small variations in performance might be observed, particularly at higher pressure or flow. As per expected clinical practice, the effectiveness of ventilation and alarms should be verified by the clinician.

⁴ Το SmartStart ενδέχεται να μην λειτουργεί σύμφωνα με την περιγραφή όταν αυτή η μάσκα χρησιμοποιείται με ορισμένες συσκευές CPAP ή δύο επιπέδων.

X/E = Η μάσκα δεν έχει εγκριθεί για τη συγκεκριμένη συσκευή.

✓ = Συμβατή X = Μη συμβατή

Μάσκες για ολόκληρο το πρόσωπο

Συσκευή	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ FX / Quattro™ FX for Her / Quattro™ Air / Quattro™ Air for Her	Mirage Liberty™	Hospital Full Face
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 11 ASV PaceWave	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AirSense™ 11	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AirMini™	X	X	X	X	X	X

Συσκευή	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ FX / Quattro™ FX for Her Quattro™ Air / Quattro™ Air for Her	Mirage Liberty™	Hospital Full Face
AirStart™ 10 / AirSense™ 10	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 CS Pacewave	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve™ 10 CS-A Pacewave	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
Lumis™	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
S9 AutoSet™ / S9 Elite™	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
S9 Escape™	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AutoSet CS™ 2	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE ¹	FULL FACE
AutoSet CS™ (S9)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AutoSet CS™ (S9 PaceWave)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
AutoSet CS-A™ (S9 PaceWave)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
C-Series Tango™		Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.				
S7™ Lightweight		Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.				
S7 Elite™ / S7 AutoSet Spirit	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S7 AutoSet Respond™	ULTRA	MIRAGE	MIRAGE	MIRAGE	MIR FULL	MIR FULL
S8™ Auto 25	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Spirit™	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Vantage™	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 Escape™		Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.				
S8 Escape™ II		Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.				
S8 Elite™ / S8 Elite™ II	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Spirit™ II	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 Compact™		Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.				
VPAP™ III / VPAP™ III ST	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ III ST-A	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ IV / VPAP™ IV ST	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ Adapt SV	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE ¹	FULL FACE
VPAP™ Adapt (S9)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Adapt (S9 PaceWave)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Malibu	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
VPAP™ Auto 25	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
Auto 25 (S9)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ S / VPAP™ ST	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ S (S9) / VPAP™ ST (S9)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ ST-A (S9)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
S9 VPAP™ Tx	X/E	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Tx	X/E	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE

Συσκευή	Mirage Series 2	Ultra Mirage™	Mirage Quattro™	Quattro™ FX / Quattro™ FX for Her Quattro™ Air / Quattro™ Air for Her	Mirage Liberty™	Hospital Full Face
Stellar™	X/E	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
VS Ultra™ / VS Integra™	X/E	✓	✓	X/E	✓	✓
VS III™	X/E	✓	✓	X/E	✓	✓
Astral™	X/E	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE

Ρινικές μάσκες 1

Συσκευή	AirTouch™ N30i ¹	AirFit™ N30 for AirMini	AirFit™ N30i ¹ / AirFit™ N30	AirFit™ N20 Classic	AirFit™ N20 / AirFit™ N20 for Her	AirTouch™ N20 / AirTouch™ N20 for Her	AirFit™ N10 / AirFit™ N10 for Her	Swift™ FX Nano / Swift™ FX Nano for Her	Mirage™ FX / Mirage FX™ for Her	Mirage Activa™ LT	Mirage™ SoftGel	Hospital Nasal
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	X	Nasal
AirCurve™ 11 ASV PaceWave	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	X	Nasal
AirSense™ 11	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	N/V	Nasal
AirMini™	X	✓	X	✓	✓	✓	X	X	X	X	X	X
AirStart™ 10 / AirSense™ 10	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 CS Pacewave	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AirCurve™ 10 CS-A Pacewave	X/E	X	X/E	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
Lumis™	Pillows ³	X	Pillows ³	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
S9 AutoSet™ / S9 Elite™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
S9 Escape™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ 2	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	VISTA ¹	ACTIVA	ACTIVA	VISTA
AutoSet CS™ (S9)	X/E	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ (S9 PaceWave)	X/E	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
AutoSet CS-A™ (S9 PaceWave)	X/E	X	X/E	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
C-Series Tango™	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.			
S7™ Lightweight	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.			
S7 Elite™ / S7 AutoSet Spirit	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	STANDARD	STANDARD	STANDARD	STANDARD
S7 AutoSet Respond™	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	STANDARD	STANDARD	STANDARD	STANDARD
S8™ Auto 25	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 AutoSet Spirit™	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 AutoSet Vantage™	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 Escape™	X/E	X	X/E	X/E	X/E	X/E	Απαιτείται βαθμολόγηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.					
S8 Escape™ II	X/E	X	X/E	X/E	Απαιτείται βαθμολόγηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.							
S8 Elite™ / S8 Elite™ II	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 AutoSet Spirit™ II	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
S8 Compact™	X/E	X	X/E	X/E	X/E	X/E	Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.					
VPAP™ III / VPAP™ III ST	X/E	X	X/E	X/E	X/E	X/E	MIRAGE	MIRAGE	STANDARD	STANDARD	STANDARD	X/E

Συσκευή	AirTouch™ N30i ¹	AirFit™ N30 for AirMini	AirFit™ N30i ¹ / AirFit™ N30	AirFit™ N20 Classic	AirFit™ N20 / AirFit™ N20 for Her	AirTouch™ N20 / AirTouch™ N20 for Her	AirFit™ N10 / AirFit™ N10 for Her	Swift™ FX Nano / Swift™ FX Nano for Her	Mirage™ FX / Mirage FX™ for Her	Mirage Activa™ LT	Mirage™ SoftGel	Hospital Nasal
VPAP™ III ST-A	X/E	X	X/E	X/E	X/E	X/E	MIRAGE	MIRAGE	STANDARD	STANDARD	STANDARD	X/E
VPAP™ IV / VPAP™ IV ST	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
VPAP™ Adapt SV	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	VISTA ¹	ACTIVA	ACTIVA	VISTA
VPAP™ Adapt (S9)	X/E	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ Adapt (S9 PaceWave)	X/E	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ Malibu	X/E	X	X/E	X/E	X/E	X/E	SWIFT	SWIFT	VISTA ¹	X/E	X/E	MERIDIAN
VPAP™ Auto 25	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
Auto 25 (S9)	X/E	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ S / VPAP™ ST	X/E	X	X/E	X/E	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	ACTIVA	ACTIVA	STANDARD
VPAP™ S (S9) / VPAP™ ST (S9)	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ ST-A (S9)	X/E	X	X/E	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	X/E	X	X/E	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
S9 VPAP™ Tx	Pillows ²	X	Pillows ²	X/E	Pillows	X/E	Pillows	Pillows	Nasal	Nasal	Nasal	Nasal
VPAP™ Tx	X/E	X	X/E	X/E	X/E	X/E	PILLOWS	PILLOWS	NASAL	NASAL	NASAL	NASAL
Stellar™	X/E	X	X/E	PILLOWS	PILLOWS	PILLOWS	PILLOWS ¹	PILLOWS ¹	NASAL	NASAL	NASAL	NASAL
VS Ultra™ / VS Integra™	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	X/E	X/E	X/E	X/E
VS III™	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E	X/E	✓	X/E	✓
Astral™	X/E	X	X/E	Pillows	Pillows	Pillows	X/E	X/E	Nasal	X/E	Nasal	Nasal

¹ Το SmartStop ενδέχεται να μην λειτουργεί αποτελεσματικά όταν αυτή η μάσκα χρησιμοποιείται με ορισμένες συσκευές CPAP ή δύο επιπέδων.

² Εξαιρούνται οι καταστάσεις λειτουργίας PAC και iVAPS. Συμβατό μόνο για το εύρος πίεσης 4 - 25 cm H2O.

³ This combination of mask and device has been determined safe for use based on the evaluation of accuracy of monitoring of pressure, volume, minute ventilation and respiratory rate, and its alarms if applicable. Small variations in performance might be observed, particularly at higher pressure or flow. As per expected clinical practice, the effectiveness of ventilation and alarms should be verified by the clinician.

X/E = Η μάσκα δεν έχει εγκριθεί για τη συγκεκριμένη συσκευή.

✓ = Συμβατή X = Μη συμβατή

Ρινικές μάσκες 2

Συσκευή	Pixi™	Mirage Kidsta™	Mirage Micro™	Mirage Vista™	Ultra Mirage II	Mirage	Mirage Activa™	Modular	Papillon™	Meridian
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	Pediatric	X	X	X	X	N/V	X	N/V	N/V	X
AirCurve™ 11 ASV PaceWave	N/V	X	X	X	X	N/V	X	N/V	N/V	X
AirSense™ 11	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V
AirMini™	X	X	X	X	X	X	X	X	X	X
AirStart™ 10 / AirSense™ 10	X/E	Nasal	Nasal	Nasal	Nasal	X/E	Nasal	X/E	X/E	Nasal
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	X/E	Nasal	Nasal	Nasal	Nasal	X/E	Nasal	X/E	X/E	Nasal
AirCurve™ 10 CS Pacewave	X/E	Nasal	Nasal	Nasal	Nasal	X/E	Nasal	X/E	X/E	Nasal
AirCurve™ 10 CS-A Pacewave	X/E	Nasal	Nasal	Nasal	Nasal	X/E	Nasal	X/E	X/E	Nasal
Lumis™	X/E	Nasal	Nasal	Nasal	Nasal	X/E	Nasal	X/E	X/E	Nasal

Συσκευή	Pixi™	Mirage Kidsta™	Mirage Micro™	Mirage Vista™	Ultra Mirage II	Mirage	Mirage Activa™	Modular	Papillon™	Meridian
S9 AutoSet™ / S9 Elite™	X/E	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
S9 Escape™	X/E	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ 2	X/E	X/E	VISTA ¹	VISTA	ULTRA	X/E	ACTIVA	X/E	PAPILLON	VISTA
AutoSet CS™ (S9)	X/E	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
AutoSet CS™ (S9 PaceWave) / AutoSet CS-A™ (S9 PaceWave)	X/E	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
C-Series Tango™	X/E	X/E	X/E				Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.			X/E
S7™ Lightweight	X/E	X/E	X/E				Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.			X/E
S7 Elite™ / S7 AutoSet Spirit	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	STANDARD	X/E
S7 AutoSet Respond™	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	STANDARD	X/E
S8™ Auto 25	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	X/E	X/E	STANDARD
S8 AutoSet Spirit™	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 AutoSet Vantage™	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 Escape™	X/E	X/E					Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.			
S8 Escape™ II	X/E	X/E					Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.			
S8 Elite™ / S8 Elite™ II	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 AutoSet Spirit™ II	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	STANDARD	STANDARD	STANDARD
S8 Compact™	X/E	X/E					Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.			
VPAP™ III / VPAP™ III ST	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	X/E	X/E
VPAP™ III ST-A	X/E	STANDARD	STANDARD	STANDARD	ULTRA	MIRAGE	STANDARD	STANDARD	X/E	X/E
VPAP™ IV / VPAP™ IV ST	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	X/E	X/E	STANDARD
VPAP™ Adapt SV	X/E	X/E	VISTA ¹	VISTA	ULTRA	X/E	ACTIVA	X/E	X/E	VISTA
VPAP™ Adapt (S9)	X/E	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
VPAP™ Adapt (S9 PaceWave)	X/E	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
VPAP™ Malibu	X/E	X/E	MERIDIAN	VISTA	ULTRA	X/E	X/E	X/E	X/E	MERIDIAN
VPAP™ Auto 25	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	X/E	X/E	STANDARD
Auto 25 (S9)	X/E	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	STANDARD
VPAP™ S / VPAP™ ST	X/E	X/E	STANDARD	STANDARD	ULTRA	MIRAGE	ACTIVA	X/E	X/E	STANDARD
VPAP™ S (S9) / VPAP™ ST (S9) / VPAP™ ST-A (S9)	Pediatric	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	Pediatric	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	Nasal	Nasal	Nasal
S9 VPAP™ Tx	Pediatric	X/E	Nasal	Nasal	Nasal Ultra	X/E	Nasal	X/E	X/E	Nasal
VPAP™ Tx	X/E	NASAL	NASAL	NASAL	NASAL-ULTRA	X/E	NASAL	X/E	X/E	NASAL
Stellar™	Pediatric	NASAL	NASAL	NASAL	ULTRA	X/E	NASAL	X/E	X/E	X/E
VS Ultra™ / VS Integra™	X/E	X/E	X/E	X/E	X/E	X/E	X/E	X/E	X/E	X/E
VS III™	X/E	✓	✓	✓	✓	X/E	✓	X/E	X/E	X/E
Astral™	Pediatric	Nasal	Nasal	Nasal	Nasal Ultra	X/E	X/E	X/E	X/E	X/E

X/E = Η μάσκα δεν έχει εγκριθεί για τη συγκεκριμένη συσκευή.

✓ = Συμβατή X = Μη συμβατή

Ρινικά μαξιλαράκια

Συσκευή	AirFit™ P30i ¹	AirFit™ P10 For AirMini™	AirFit™ P10 / AirFit™ P10 for Her	Mirage Swift™	Mirage Swift™ II	Swift™ FX / Swift™ FX for Her / Swift™ FX Bella / Swift™ FX Bella Gray	Swift™ LT	Swift™ LT for Her / Swift™ LT-F
AirCurve™ 11 VAuto / AirCurve™ 11 S / AirCurve™ 11 ST	Pillows	X	Pillows ²	X/E	X/E	Pillows	Pillows	Pillows
AirCurve™ 11 ASV PaceWave	Pillows	X	Pillows ²	X/E	X/E	Pillows	Pillows	Pillows
AirSense™ 11	Pillows	X	Pillows	X/E	Pillows	Pillows	Pillows	Pillows
AirMini™	X	✓	X	X	X	X	X	X
AirStart™ 10 / AirSense™ 10	Pillows	X	Pillows ²	X/E	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	Pillows	X	Pillows ²	X/E	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 CS Pacewave	Pillows	X	Pillows ²	X/E	Pillows	Pillows	Pillows	Pillows
AirCurve™ 10 CS-A Pacewave	X/E	X	Pillows ²	X/E	Pillows	Pillows	Pillows	Pillows
Lumis™	Pillows ⁴	X	Pillows ²	X/E	Pillows	Pillows	Pillows	Pillows
S9 AutoSet™ / S9 Elite™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows
S9 Escape™	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows	Pillows
AutoSet CS™ 2	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E
AutoSet CS™ (S9)	Pillows	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
AutoSet CS™ (S9 PaceWave)	Pillows	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
AutoSet CS-A™ (S9 PaceWave)	X/E	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
C-Series Tango™	X/E	X				Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.		
S7™ Lightweight	X/E	X				Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.		
S7 Elite™ / S7 AutoSet Spirit	X/E	X	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
S7 AutoSet Respond™	X/E	X	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
S8™ Auto 25	X/E	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S8 AutoSet Spirit™	X/E	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S8 AutoSet Vantage™	X/E	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S8 Escape™	X/E	X				Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.		
S8 Escape™ II	X/E	X				Απαιτείται βαθμονόμηση της πίεσης στα 12 cm H2O. Συμβουλευτείτε τον ιατρό σας και ελέγξτε τις οδηγίες χρήσης.		
S8 Elite™ / S8 Elite™ II	X/E	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S8 AutoSet Spirit™ II	X/E	X	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
S8 Compact™	X/E	X				Δεν απαιτείται ρύθμιση. Συμβουλευτείτε τον ιατρό σας πριν χρησιμοποιήσετε διαφορετικό τύπο μάσκας.		
VPAP™ III / VPAP™ III ST	X/E	X	MIRAGE ²	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
VPAP™ III ST-A	X/E	X	MIRAGE ²	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
VPAP™ IV / VPAP™ IV ST	X/E	X	SWIFT ²	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ Adapt SV	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E
VPAP™ Adapt (S9)	Pillows	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ Adapt (S9 PaceWave)	Pillows	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ Malibu	X/E	X	SWIFT ²	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ Auto 25	X/E	X	SWIFT ²	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
Auto 25 (S9)	Pillows	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows

Συσκευή	AirFit™ P30i ¹	AirFit™ P10 For AirMini™	AirFit™ P10 / AirFit™ P10 for Her	Mirage Swift™	Mirage Swift™ II	Swift™ FX / Swift™ FX for Her / Swift™ FX Bella / Swift™ FX Bella Gray	Swift™ LT	Swift™ LT for Her / Swift™ LT-F
VPAP™ S / VPAP™ ST	X/E	X	SWIFT ²	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT
VPAP™ S (S9) / VPAP™ ST (S9)	Pillows	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ ST-A (S9)	X/E	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
VPAP™ ST (S9 iVAPS) / VPAP™ ST-A (S9 iVAPS)	X/E	X	Pillows ²	Pillows	Pillows	Pillows	Pillows	Pillows
S9 VPAP™ Tx	Pillows ³	X	Pillows ²	X/E	Pillows	Pillows	Pillows	Pillows
VPAP™ Tx	X/E	X	PILLOWS ²	PILLOWS	PILLOWS	PILLOWS	PILLOWS	PILLOWS
Stellar™	X/E	X	PILLOWS ^{1,2}	PILLOWS ¹	PILLOWS ¹	PILLOWS ¹	PILLOWS ¹	PILLOWS ¹
VS Ultra™ / VS Integra™	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E
VS III™	X/E	X	X/E	X/E	X/E	X/E	X/E	X/E
Astral™	X/E	X	X/E	Pillows	Pillows	Pillows	Pillows	Pillows

¹ Το SmartStop ενδέχεται να μην λειτουργεί αποτελεσματικά όταν αυτή η μάσκα χρησιμοποιείται με ορισμένες συσκευές CPAP ή δύο επιπέδων.

² Η ακρίβεια του αναπνεόμενου όγκου/αερισμού κατά λεπτό ενδέχεται να μειωθεί κατά τη χρήση της μάσκας AirFit P10 σε λειτουργίες δύο επιπέδων (bilevel) ειδικά στις υψηλότερες υποστηρίξεις πίεσης.

³ Εξαιρούνται οι καταστάσεις λειτουργίας PAC και iVAPS. Συμβατό μόνο για το εύρος πίεσης 4 - 25 cm H₂O.

⁴ This combination of mask and device has been determined safe for use based on the evaluation of accuracy of monitoring of pressure, volume, minute ventilation and respiratory rate, and its alarms if applicable. Small variations in performance might be observed, particularly at higher pressure or flow. As per expected clinical practice, the effectiveness of ventilation and alarms should be verified by the clinician.

X/E = Η μάσκα δεν έχει εγκριθεί για τη συγκεκριμένη συσκευή.

✓ = Συμβατή X = Μη συμβατή